

CUSTOMER SERVICE BULLETIN

Radio Shack® CUSTOMER SERVICE BULLETIN

400 Atrium—One Tandy Center, Fort Worth, Texas 76102

Bulletin No.	1568-Ø7
Description	MEDICAL OFFICE
Model	Model I/III
Stock No.	26-1568 Version 1.Ø
Date	August 9, 1982
Page 1 of	1

Problem

Error 15 in Line 113Ø of Patient Sort/Print: This problem will occur only when the combination of a title (i.e. Miss, Rev, Dr) is entered as the first name, as in 'Jones, Dr' with no characters after the title, and a space between the comma and the title. Other combinations will not cause the problem, as in 'Jones,Dr' or 'Jones, Miss E.'

* * * THESE LINE CHANGES ARE OPTIONAL! * * *

NOTE: These line changes are optional, as they will slow the sort down in checking for an extremely unusual situation.

Solution

The problem is corrected by following the steps listed below.

1. Backup the diskette(s) and make the changes on the Backup copy of the program.
2. In BASIC load the program by typing LOAD"ALPHASOR".
3. Make the following corrections:

CHANGES (Retype the line or refer to the Edit section of the owners manual)

Old Line: 113Ø IFINSTR(MID\$(NV\$,J+1),PR\$(II))=1THENNV\$=MID\$(NV\$,1,J)+
MID\$(NV\$,J+1+LEN(PR\$(II)))+ " "+PR\$(II):II=Ø:GOTO113ØELSENEXTII

New Line: 113Ø IFINSTR(MID\$(NV\$,J+1),PR\$(II))=1THENNV\$=MID\$(NV\$,1,J)+
MID\$(NV\$,J+1+LEN(PR\$(II)))ELSENEXTII:GOTO114Ø

ADD THIS LINE: 1135 IFII=ØTHEN113ØELSENV\$=NV\$+" "+PR\$(II):II=Ø:GOTO113Ø

4. Type SAVE"ALPHASOR" to save the changes in the program.
5. At TRSDOS READY, make a backup copy of the corrected diskette.